

Verifying Readiness of the CR2

If the device is not connected to the WIFI network you should check the green readiness indicator on the device at least once each month or as per local protocol. If the device is not ready, the readiness indicator does not flash and an alert tone sounds every 15 minutes.

If the device is not ready, perform the following steps to determine why the readiness indicator is not flashing.

1. Open the lid and **wait** until the voice prompts start.



2. **Immediately** press and hold the **LANGUAGE** and **CHILD MODE** buttons simultaneously until you hear either **DEVICE READY** or **DEVICE NOT READY**.



3. The defibrillator will then provide voice instructions that tell you which of the following actions you need to perform:
 - Replace electrode tray
 - Replace battery
 - Contact qualified service personnel

Note: If you pressed both buttons but did not hear the **DEVICE READY** or **DEVICE NOT READY** voice prompts, you may need to try again. When you open the lid, you must wait until the voice prompts start before pressing the two buttons. After the voice prompts start, you must press the two buttons within 10 seconds. If you miss these time limits, the AED proceeds as it would during a cardiac arrest. To try again, close and reopen the lid.