

Welcome to PowerMic Mobile

PowerMic Mobile enables you to use your mobile device as a microphone for Dragon Medical One Desktop Application. PowerMic Mobile on your device pairs with Dragon Medical One on your PC.



Start Dragon Medical One, enter your User name and select Nuance **PowerMic Mobile** from the **Microphone list**.



You will be presented with the Dragon Bar 'Amber' waiting for connection Icon.



You will firstly need to download Powermic Mobile from the App store / Play store.



Your mobile device and your PC must have a network connection so they can communicate with each other. Proceed as follows:

1. Copy the configuration link (ensure you have the correct link)

Config links


-----**Android config link**-----

http://config_?NmsToken=MzU5MWFIM2ItYjUwOS00OWUzLTk4YmEtYmU3MDA4MjNlMjVI&NmsBaseUrl=nms-uk.nuancehdp.com

-----**iOS config link**-----

dmic://config_?NmsToken=MzU5MWFIM2ItYjUwOS00OWUzLTk4YmEtYmU3MDA4MjNlMjVI&NmsBaseUrl=nms-uk.nuancehdp.com



2. Start PowerMic Mobile and tap Add Profile



Click yes to paste the URL as per the above pic. You can also do this manually by copying the URL link into the PROFILE URL section.


Remember to use the correct config string that relates to the type of smart phone you are using (iOS or Android)

You will notice that once pasted, 'Valid profile' will be displayed as per the below image.



Enter a profile Name. This does NOT have to be your User ID, but can also be a friendly name, like 'Andrews PowerMic'



Click 'Save'



Now Enter your User name.


This will be the same name you logged into your Dragon Medical One Client.

Click 'Log In'




After you get the 'Welcome' prompt, click 'Start', then Log In.

Once successfully connected, your controls will be displayed as per below:



Your Dragon Medical Client 'Dragon bar' Icon will also turn Red.

You are ready now to dictate. Click on the Blue icon to record. Your Icon will turn a solid blue colour. Press again to pause. Your icon will invert.




Your Dragon Medical One Icon will also turn green to indicate you are in record.

Training a word using the PowerMic Mobile.


Training a particular word in PowerMic Mobile can be somewhat challenging.

Try using the 'Enter the word as you would pronounce it' whenever possible.


1. Enter the word you wish to train. (The example shows adding the word 'PowerMic Mobile')
2. Click on the Blue record button on your iPhone/Android (it should go solid blue)
3. Click on the Red icon on your Dragon Medical Client.



4. The icon will turn green and start revolving. At this point, say the word you added. Say this in the most natural voice without enunciating the word.




5. Once your word has been recognised, you can click the 'Confirm' icon.
6. Now turn your Microphone off, by clicking the Blue button on your iPhone/Android device.



7. Your success will be confirmed. (see below)

Add Word



**If the word was trained in the incorrect sequence, the following error could occur. If this happens, try again from the beginning, by clicking 'retry'.

